



NASO-FLO[®]

*Naso-Flo[®] Nasopharyngeal
Airway Device*



The Naso-Flo[®] nasopharyngeal airway device allows for direct oxygen delivery, while humidification vents positioned towards the distal tip facilitate heat and moisture transfer. It can also be supplied with an optional respiratory indicator with or without hydrophobic filter.

- Manufactured from soft silicone **reducing the risk of trauma** during insertion
- **Colour-coded connectors** to aid correct size selection
- **Delivers oxygen directly into the pharynx**, reducing the anatomical dead space
- Once the patient's tidal volume has been exceeded, the **upper airway becomes a natural reservoir** for the remaining oxygen
- Available with or without a Respiratory Indicator, which **allows for detection of CO₂ from the distal tip**
- **Optional hydrophobic filter** (on versions with Respiratory Indicator), which reduces the risk of monitor contamination
- Suitable for a **wide range of clinical applications**
- **Humidification vents** and the physiological position improves the humidification of oxygen
- **Ventilation adaptor** available for connection between Naso-Flo[®] and breathing system
- **Cost effective** - Reduces the need for an oxygen mask with the additional benefit of improving the visibility of the patients oral cavity
- **MRI Compatible**
- **Single Use** 
- **Latex Free** 

FOR MORE INFORMATION ON OUR PRODUCTS AND SERVICES PHONE US ON +44 (0)7908 789 712

Medis Medical Tianjin Company Ltd specialises in the design, development, sales and marketing of unique anaesthesia and airway management devices worldwide. Led by a team of pioneering inventors, we have many years experience bringing new products to market that solve problems, improve patient safety and change lives.



ORDERING INFORMATION

With O2 Port

Size	Code	Box Qty
4.0mm	NF104	10
5.0mm	NF105	10
6.0mm	NF106	10
7.0mm	NF107	10
8.0mm	NF108	10
9.0mm	NF109	10



With O2 Port, Respiratory Indicator and Filter

Size	Code	Box Qty
4.0mm	NF204	10
5.0mm	NF205	10
6.0mm	NF206	10
7.0mm	NF207	10
8.0mm	NF208	10
9.0mm	NF209	10



With Respiratory Indicator, without Filter

Size	Code	Box Qty
4.0mm	NF304	10
5.0mm	NF305	10
6.0mm	NF306	10
7.0mm	NF307	10
8.0mm	NF308	10
9.0mm	NF309	10



Ventilation adaptor

Size	Code	Box Qty
15mm	NF400	50

