



OXYCAP™

OxyCap™ universal recovery adaptor



The OxyCap™ universal recovery adaptor allows for direct oxygen delivery into the airway tube. It is suitable for use with all types of Laryngeal Masks and Tracheostomy devices. Available in two versions, with and without respiratory indicator.

- **Available in two versions** – with and without respiratory indicator. With respiratory indicator allows for detection of CO₂ from within the airway tube
- Cost effective – **reduces the need for an oxygen mask** by facilitating direct oxygen delivery into the airway tube
- **Suitable for use on all types** of Laryngeal Masks & Tracheostomy devices on spontaneously breathing patients
- **Red colour connector** – allows for quick, visual identification of OxyCap™ connected patients
- **Hydrophobic Filter** – fitted to the respiratory indicator and reduces the risk of monitor contamination
- **Lightweight device** – reduces the risk of accidental displacement of the airway tube
- **Individually Wrapped**
- **MRI Compatible**
- **Single Use** 
- **Latex Free** 

FOR MORE INFORMATION ON OUR PRODUCTS AND SERVICES PHONE US ON +44 (0)7908 789 712

Medis Medical Tianjin Company Ltd specialises in the design, development, sales and marketing of unique anaesthesia and airway management devices worldwide. Led by a team of pioneering inventors, we have many years experience bringing new products to market that solve problems, improve patient safety and change lives.



ORDERING INFORMATION

Product	Code
OxyCap™ with respiratory indicator	OC0100
OxyCap™ without respiratory indicator	OC0200

