




PAT-TUBE®

Pat-Tube® endotracheal tube with integrated throat/pharyngeal pack



Pat-Tube® is a specialised endotracheal tube with an integrated throat/pharyngeal pack that is used instead of a loose throat pack, eliminating the risk of accidental retention. It offers clinicians a safe, clinical alternative whenever a throat pack is required and an Endotracheal Tube is used.

- High volume, low pressure pharyngeal pack **eliminates the issue of throat pack retention and reduces the risk of sore throat** following surgery
- Inflatable pharyngeal pack **delivers a complete fluid seal**, preventing surgical debris and blood entering the lower airway or oesophagus
- Oral and nasal versions allow for a **wide variety of clinical applications**
- Supplied with 2 pressure indicators **enabling continuous visual monitoring** of both the pharyngeal pack and tracheal cuff pressure during use
- Oral versions are supplied with a reinforced tube, **reducing the risk of occlusion**
- Preformed north facing nasal tube allows for **greater access** during oral surgical procedures
- Also available in **reinforced and preformed** nasal & oral versions
- **Paediatric sizes** available in most versions
- All versions supplied with **Murphy Eye**
- **Single Use** 
- **Latex Free** 
- **DEHP Free** 

FOR MORE INFORMATION ON OUR PRODUCTS AND SERVICES PHONE US ON +44 (0)7908 789 712

Medis Medical Tianjin Company Ltd specialises in the design, development, sales and marketing of unique anaesthesia and airway management devices worldwide. Led by a team of pioneering inventors, we have many years experience bringing new products to market that solve problems, improve patient safety and change lives.



ORDERING INFORMATION

ADULT

Pat-Tube®, Standard, Oral/Nasal		
Size	Code	Box Qty
6.0mm	60060	10
6.5mm	60065	10
7.0mm	60070	10
7.5mm	60075	10
8.0mm	60080	10
8.5mm	60085	10
9.0mm	60090	10

Pat-Tube®, Preformed, South Facing, Oral		
Size	Code	Box Qty
6.0mm	50060	10
6.5mm	50065	10
7.0mm	50070	10
7.5mm	50075	10
8.0mm	50080	10
8.5mm	50085	10
9.0mm	50090	10

Pat-Tube®, Preformed, North Facing, Nasal		
Size	Code	Box Qty
6.0mm	20060	10
6.5mm	20065	10
7.0mm	20070	10
7.5mm	20075	10

Pat-Tube®, Reinforced, Oral		
Size	Code	Box Qty
6.0mm	10060	10
6.5mm	10065	10
7.0mm	10070	10
7.5mm	10075	10
8.0mm	10080	10
8.5mm	10085	10
9.0mm	10090	10

PAEDIATRIC

Pat-Tube®, Reinforced, Oral, Uncuffed		
Size	Code	Box Qty
4.5mm	30045	10
5.0mm	30050	10
5.5mm	30055	10

Pat-Tube®, Reinforced, Oral, Uncuffed		
Size	Code	Box Qty
4.5mm	40045	10
5.0mm	40050	10
5.5mm	40055	10

