





TRACHEOSTOMY TUBES

*Adult and paediatric
tracheostomy tubes*



Tracheostomy Tubes are supplied cuffed and uncuffed with optional AccuCuff™ cuff pressure indicator to suit a variety of clinical situations and patient needs.

- **Standard and reinforced** versions available
- Comprehensive range of **cuffed and uncuffed tubes**
- Optional AccuCuff™ cuff pressure indicator allows for continuous visual monitoring, eliminating the need for regular spot-checking and **reducing the risk of Post Ventilatory Pneumonia**
- Negative pressure graduation on the AccuCuff™ cuff pressure indicator ensures cuff is fully deflated, **reducing the risk of trauma during removal**
- Fully integrated inflation lines connect to the Tracheostomy Tube proximal to the patient connector **reducing the risk of chafing and infection to the Tracheostomy site**
- **Single Use** 
- **Latex Free** 
- **DEHP Free** 
- Supplied in **boxes of 10**

FOR MORE INFORMATION ON OUR PRODUCTS AND SERVICES PHONE US ON +44 (0)7908 789 712

Medis Medical Tianjin Company Ltd specialises in the design, development, sales and marketing of unique anaesthesia and airway management devices worldwide. Led by a team of pioneering inventors, we have many years experience bringing new products to market that solve problems, improve patient safety and change lives.



ORDERING INFORMATION

Tracheostomy Tube

Tracheostomy Tube (uncuffed)				Tracheostomy Tube with cuff				Tracheostomy Tube with AccuCuff™			
Size	Code	Size	Code	Size	Code	Size	Code	Size	Code	Size	Code
2.0mm	TT-P1-20	6.5mm	TT-P1-65	2.0mm	-	6.5mm	TT-P2-65	2.0mm	-	6.5mm	TT-P3-65
2.5mm	TT-P1-25	7.0mm	TT-P1-70	2.5mm	-	7.0mm	TT-P2-70	2.5mm	-	7.0mm	TT-P3-70
3.0mm	TT-P1-30	7.5mm	TT-P1-75	3.0mm	TT-P2-30	7.5mm	TT-P2-75	3.0mm	TT-P3-30	7.5mm	TT-P3-75
3.5mm	TT-P1-35	8.0mm	TT-P1-80	3.5mm	TT-P2-35	8.0mm	TT-P2-80	3.5mm	TT-P3-35	8.0mm	TT-P3-80
4.0mm	TT-P1-40	8.5mm	TT-P1-85	4.0mm	TT-P2-40	8.5mm	TT-P2-85	4.0mm	TT-P3-40	8.5mm	TT-P3-85
4.5mm	TT-P1-45	9.0mm	TT-P1-90	4.5mm	TT-P2-45	9.0mm	TT-P2-90	4.5mm	TT-P3-45	9.0mm	TT-P3-90
5.0mm	TT-P1-50	9.5mm	TT-P1-95	5.0mm	TT-P2-50	9.5mm	TT-P2-95	5.0mm	TT-P3-50	9.5mm	TT-P3-95
5.5mm	TT-P1-55	10.0mm	TT-P1-100	5.5mm	TT-P2-55	10.0mm	TT-P2-100	5.5mm	TT-P3-55	10.0mm	TT-P3-100
6.0mm	TT-P1-60			6.0mm	TT-P2-60			6.0mm	TT-P3-60		

Tracheostomy Tube Reinforced

Tracheostomy Tube Reinforced (uncuffed)				Tracheostomy Tube Reinforced with cuff				Tracheostomy Tube Reinforced with AccuCuff™			
Size	Code	Size	Code	Size	Code	Size	Code	Size	Code	Size	Code
2.0mm	-	6.5mm	TT-J1-65	2.0mm	-	6.5mm	TT-J2-65	2.0mm	-	6.5mm	TT-J3-65
2.5mm	-	7.0mm	TT-J1-70	2.5mm	-	7.0mm	TT-J2-70	2.5mm	-	7.0mm	TT-J3-70
3.0mm	TT-J1-30	7.5mm	TT-J1-75	3.0mm	TT-J2-30	7.5mm	TT-J2-75	3.0mm	TT-J3-30	7.5mm	TT-J3-75
3.5mm	TT-J1-35	8.0mm	TT-J1-80	3.5mm	TT-J2-35	8.0mm	TT-J2-80	3.5mm	TT-J3-35	8.0mm	TT-J3-80
4.0mm	TT-J1-40	8.5mm	TT-J1-85	4.0mm	TT-J2-40	8.5mm	TT-J2-85	4.0mm	TT-J3-40	8.5mm	TT-J3-85
4.5mm	TT-J1-45	9.0mm	TT-J1-90	4.5mm	TT-J2-45	9.0mm	TT-J2-90	4.5mm	TT-J3-45	9.0mm	TT-J3-90
5.0mm	TT-J1-50	9.5mm	TT-J1-95	5.0mm	TT-J2-50	9.5mm	TT-J2-95	5.0mm	TT-J3-50	9.5mm	TT-J3-95
5.5mm	TT-J1-55	10.0mm	TT-J1-100	5.5mm	TT-J2-55	10.0mm	TT-J2-100	5.5mm	TT-J3-55	10.0mm	TT-J3-100
6.0mm	TT-J1-60			6.0mm	TT-J2-60			6.0mm	TT-J3-60		

