





# CUFF-SAFE®

*Cuff-Safe® endotracheal tube*



**Cuff-Safe® endotracheal tube range, supplied with AccuCuff™ cuff pressure indicator as standard.**

- **Supplied with** AccuCuff™ cuff pressure indicator
- AccuCuff™ cuff pressure indicator allows for continuous visual monitoring of the cuff pressure, **eliminating the need for regular spot-checking**
- Negative pressure graduation on the AccuCuff™ indicates when the cuff is fully deflated, **reducing the risk of trauma during removal**
- **Extensive range** of standard, reinforced, preformed cuffed tubes available
- Manufactured from **transparent pvc**
- Reinforced versions available **reducing the risk of accidental occlusion**
- **Clear connectors – EGBAT compliant**
- **Murphy Eye**
- **Single Use** 
- **Latex Free** 
- **DEHP Free** 
- 5-year shelf life allows **long term storage and reduces waste**
- Supplied in **boxes of 10**

**FOR MORE INFORMATION ON OUR PRODUCTS AND SERVICES PHONE US ON +44 (0)7908 789 712**

Medis Medical Tianjin Company Ltd specialises in the design, development, sales and marketing of unique anaesthesia and airway management devices worldwide. Led by a team of pioneering inventors, we have many years experience bringing new products to market that solve problems, improve patient safety and change lives.



# ORDERING INFORMATION

## Standard Endotracheal Tube

Standard Endotracheal Tube with AccuCuff™			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-P31-65
2.5mm	ETT-P31-25	7.0mm	ETT-P31-70
3.0mm	ETT-P31-30	7.5mm	ETT-P31-75
3.5mm	ETT-P31-35	8.0mm	ETT-P31-80
4.0mm	ETT-P31-40	8.5mm	ETT-P31-85
4.5mm	ETT-P31-45	9.0mm	ETT-P31-90
5.0mm	ETT-P31-50	9.5mm	ETT-P31-95
5.5mm	ETT-P31-55	10.0mm	ETT-P31-100
6.0mm	ETT-P31-60		

Standard Endotracheal Tube with AccuCuff™ (High Volume, Low Pressure Cuff)			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-P32-65
2.5mm	ETT-P32-25	7.0mm	ETT-P32-70
3.0mm	ETT-P32-30	7.5mm	ETT-P32-75
3.5mm	ETT-P32-35	8.0mm	ETT-P32-80
4.0mm	ETT-P32-40	8.5mm	ETT-P32-85
4.5mm	ETT-P32-45	9.0mm	ETT-P32-90
5.0mm	ETT-P32-50	9.5mm	ETT-P32-95
5.5mm	ETT-P32-55	10.0mm	ETT-P32-100
6.0mm	ETT-P32-60		

## Endotracheal Tube Oral Preformed South Facing

Endotracheal Tube Oral Preformed with AccuCuff™			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-K31-65
2.5mm	-	7.0mm	ETT-K31-70
3.0mm	ETT-K31-30	7.5mm	ETT-K31-75
3.5mm	ETT-K31-35	8.0mm	ETT-K31-80
4.0mm	ETT-K31-40	8.5mm	ETT-K31-85
4.5mm	ETT-K31-45	9.0mm	ETT-K31-90
5.0mm	ETT-K31-50	9.5mm	ETT-K31-95
5.5mm	ETT-K31-55	10.0mm	ETT-K31-100
6.0mm	ETT-K31-60		

Endotracheal Tube Oral Preformed with AccuCuff™ (High Volume, Low Pressure Cuff)			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-K32-65
2.5mm	-	7.0mm	ETT-K32-70
3.0mm	ETT-K32-30	7.5mm	ETT-K32-75
3.5mm	ETT-K32-35	8.0mm	ETT-K32-80
4.0mm	ETT-K32-40	8.5mm	ETT-K32-85
4.5mm	ETT-K32-45	9.0mm	ETT-K32-90
5.0mm	ETT-K32-50	9.5mm	ETT-K32-95
5.5mm	ETT-K32-55	10.0mm	ETT-K32-100
6.0mm	ETT-K32-60		

## Endotracheal Tube Nasal Preformed North Facing

Endotracheal Tube Nasal Preformed with AccuCuff™			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-B31-65
2.5mm	-	7.0mm	ETT-B31-70
3.0mm	ETT-B31-30	7.5mm	ETT-B31-75
3.5mm	ETT-B31-35	8.0mm	ETT-B31-80
4.0mm	ETT-B31-40	8.5mm	ETT-B31-85
4.5mm	ETT-B31-45	9.0mm	ETT-B31-90
5.0mm	ETT-B31-50	9.5mm	ETT-B31-95
5.5mm	ETT-B31-55	10.0mm	ETT-B31-100
6.0mm	ETT-B31-60		

Endotracheal Tube Nasal Preformed with AccuCuff™ (High Volume, Low Pressure Cuff)			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-B32-65
2.5mm	-	7.0mm	ETT-B32-70
3.0mm	ETT-B32-30	7.5mm	ETT-B32-75
3.5mm	ETT-B32-35	8.0mm	ETT-B32-80
4.0mm	ETT-B32-40	8.5mm	ETT-B32-85
4.5mm	ETT-B32-45	9.0mm	ETT-B32-90
5.0mm	ETT-B32-50	9.5mm	ETT-B32-95
5.5mm	ETT-B32-55	10.0mm	ETT-B32-100
6.0mm	ETT-B32-60		

## Endotracheal Tube Reinforced

Endotracheal Tube Reinforced with AccuCuff™			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-J31-65
2.5mm	-	7.0mm	ETT-J31-70
3.0mm	ETT-J31-30	7.5mm	ETT-J31-75
3.5mm	ETT-J31-35	8.0mm	ETT-J31-80
4.0mm	ETT-J31-40	8.5mm	ETT-J31-85
4.5mm	ETT-J31-45	9.0mm	ETT-J31-90
5.0mm	ETT-J31-50	9.5mm	ETT-J31-95
5.5mm	ETT-J31-55	10.0mm	ETT-J31-100
6.0mm	ETT-J31-60		

Endotracheal Tube Reinforced with AccuCuff™ (High Volume, Low Pressure Cuff)			
Size	Code	Size	Code
2.0mm	-	6.5mm	ETT-J32-65
2.5mm	-	7.0mm	ETT-J32-70
3.0mm	ETT-J32-30	7.5mm	ETT-J32-75
3.5mm	ETT-J32-35	8.0mm	ETT-J32-80
4.0mm	ETT-J32-40	8.5mm	ETT-J32-85
4.5mm	ETT-J32-45	9.0mm	ETT-J32-90
5.0mm	ETT-J32-50	9.5mm	ETT-J32-95
5.5mm	ETT-J32-55	10.0mm	ETT-J32-100
6.0mm	ETT-J32-60		

