




# BOUGIE

*tracheal tube bougies*



bougies are designed to aid the placement of endotracheal and tracheostomy tubes during difficult intubations. They are supplied with a range of features that allow for easier insertion, reducing trauma and improving patient comfort.

- Internal memory allowing for **improved positioning during insertion**
- **High visibility colour** to aid positioning in poor light
- Low point of pressure tip **reducing the risk of inadvertent trauma during insertion**
- **Adult and Paediatric versions available** for use in a wide range of clinical applications
- Available in 5ch, 10ch and 15ch diameter and in 600mm or 800mm lengths
- Both **angled (coudé) tip and straight versions** available
- Low friction outer surface **enhances tracheal tube glide**
- **Single Use** 
- **Latex Free** 
- Supplied in **boxes of 10**

**FOR MORE INFORMATION ON OUR PRODUCTS AND SERVICES PHONE US ON +86-22-83963862**

Medis Medical Tianjin Company Ltd specialises in the design, development, sales and marketing of unique anaesthesia and airway management devices worldwide. Led by a team of pioneering inventors, we have many years experience bringing new products to market that solve problems, improve patient safety and change lives.



# ORDERING INFORMATION

Code	Description	Size
47005S	Bougie, straight	5ch x 470mm
60010S	Bougie, straight	10ch x 600mm
60015S	Bougie, straight	15ch x 600mm
80010S	Bougie, straight	10ch x 800mm
80015S	Bougie, straight	15ch x 800mm
60010A	Bougie, angled (coudé)	10ch x 600mm
60015A	Bougie, angled (coudé)	15ch x 600mm
80010A	Bougie, angled (coudé)	10ch x 800mm
80015A	Bougie, angled (coudé)	15ch x 800mm



Internal Memory